Basic Series

The orthopedic insoles of the Basic Series Offer great support with a plastic base and a neutral intrinsic post.

The semi-rigid shell can be reinforced using an arch fill under the plate with soft foam and keep some flexibility.

The intrinsic post reduces bulk and stabilizes the foot.





Orthopedic Insoles—Basic Series



Legionorthopedics.com

Legion Orthopedics LLP

P.O. BOX 67

Deerfield street, NJ 08313 USA

Tel:(856)3925834 Fax:(856)4595187 info@legionorthopedics.com


Authorized Distributor:



WHITE - BLACK POLYPROPYLENE

Shell - Base measurements

● ● ■ LEGION ORTHOPEDICS—Basic Series




Use this chart to know the approximate measurements of the white polypropylene shell of the Basic Series Insoles

The final size and shell election should be approved by the doctor or foot care provider

Please note this are pre made foot orthotics made to an average arch foot, may not fit everybody. Please check our different series.

Use the next picture as reference

- 1. Arch side length (medial)
- 2. Lateral length
- 3. Heel width
- 4. Distal width



● ● ● Shell - Base measurements / Sizes approximations

LETTER	1	2	3	4	
	Length*	Length*	Heel*	Distal*	
Α	15.2	13.7	6.2	7.2	
В	15.8	14.0	6.4	7.3	
С	16.5	14.9	6,6	7.5	
D	17.0	15.7	6.9	7.7	
Е	17.5	16.2	7.1	7.9	
F	18.2	16.9	7.3	8.2	

^{*}Measurements are in centimeters



Pic 1 represents recommended length (close to metatarsal heads)

Pic 2 The shell is a bit short

Pic 3 The shell is too long

Recommended Shoe Size (US)

Letter	Α	В	С	D	E	F
Man	5- 6.5	7- 7.5	8- 8.5	9- 9.5	10- 10. 5	11- 13
Woman	6- 7.5	8- 8.5	9- 9.5	10- 10. 5	11- 11. 5	12- 14

Note: this recommended sizes should be changed if necessary, adjusted accordingly to better fit the arch length of your patient's foot. Use the order form with the Basic Series Kit if available for better fitment.

Legion Orthopedics LLP

P.O. BOX 67

Deerfield street, NJ 08313 USA

Tel:(856)3925834 Fax:(856)4595187 info@legionorthopedics.com

Authorized Distributor: